permedica s.p.a. P.IVA: IT 01975020130 Via Como 38/39, 23807 Merate (Lc) - ITALY Tel.: +39 039 9514811 Fax: +39 039 9903078 E-mail: info@permedica.it web: www.permedica.it

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1. DEFINITION OPL - Osseointegration Prosthetic Limb system: set of implantable and non-implantable devices, suitable for intramedullary and percutaneous bone fixation of an armal lower or upper limb prosthesis in patients with amputation.

2. GENERAL INFORMATION Before using Permedica Spa products, the surgeon must carefully study the information below and the information available relating to the specific product (Technical Data Sheet and Surgical Technique).

Failure to comply with the instructions reported in this instruction sheet exempts the manufacturer from any liability. For clinical performances and clinical safety claims please check the summary of safety and clinical performances (SSCP) available in EUADAMED (European database on medical devices) or available upon request to Permedica Spa.

2.1 INTENDED PURPOSE, PERFORMANCES CHARACTERISTICS AND CLINICAL BENEFITS The OPL system is intended for the percutaneous invasive treatment of patients with lower or upper limb amputation. The purpose of the OPL system is to provide, through a percutaneous system, a stable and osseointegrated fixation between the bone stump of the amputated limb and an external artificial limb prosthesis.

The OPL system provides the amputee patient a stable and lasting external limb prosthesis directly attached to the remaining bone stump of the amputated limb, resulting in better proprioception and use of the artificial limb, eliminating complications related to the use of a traditional external prosthesis, thus allowing the patient to return to normal daily life activities.

The OPL system, however, even if successfully implanted, will never be able to replace a natural and healthy limb. The expected duration for the proper functioning of the OPL system is limited by residual risks associated with its use, that could result in failure and subsequent surgical removal. If the user of the OPL system and the patient comply with the precautions and warnings below, along with adequate postoperative care, it is possible to minimise the residual risks and provide a usage lifetime of many years.

2.2 PATIENT TARGET GROUP Devices of the OPL system are intended for adult patients with mature skeleton (closed physes), who have undergone amputation of the lower limb at the trans-femoral or trans-tibial level, or amputations of the upper limb at the trans-humeral level.

2.3 INTENDED USERS The intended users of the devices of the OPL system are orthopedic surgeons qualified in traumatology and joint replacement surgery. In particular, specialists in the percutaneous osseointegration of external prosthetic systems of artificial limbs.

3. MATERIALS The materials used for the production of each implantable device are shown on the label. The devices can be produced with the following materials:

Table with 2 columns: Material, Reference Standard. Rows include Titanium alloy Ti6Al7Nb (ISO 5832/11) and CrCoMo forged alloy (ISO 5832/12).

Table with 2 columns: Material, Reference Standard. Rows include Pure Titanium (ASTM F1580) and TiNiNb ( / ).

Not implantable devices are produced in titanium alloys, steel and plastic materials. 4. INFORMATION AND PRECAUTIONS FOR THE SAFE USE OF THE DEVICES 4.1 SYSTEM DESCRIPTION Permedica Spa's OPL fixation system is a modular system intended only for percutaneous intramedullary bone fixation through which an external limb prosthesis (exoprosthesi) is connected directly to the remaining bone stump in amputee patients.

The OPL fixation system consists of a series of modular components that can be assembled. The system includes both implantable and non-implantable (non-invasive) components and consists of: ➢ an intramedullary stem, which constitutes the component directly connected to the bone of the prosthetic fixation system by osteointegration. The stems are available in four different versions, differing in curvature, shapes and connection type;

- an end cap for temporary protection of the stem;
• an end cap for tightening the percutaneous coupling to the stem;
• an external coupling to attach/detach the limb-exoprosthesi.

The prosthetic components of the OPL system have been developed to give exoprosthesi of the lower or upper limb a direct anchorage with the remaining portion of the bone stump in patients with amputations and have been designed to mechanically support efforts related to a full load for a limited period of time.

In particular, the in situ lifetime of the double cone percutaneous coupling is 2 years, unless there are any previous contraindications. In general, the duration of proper functioning of the bone fixation system for exoprosthesi of limbs is limited, as it is subject to inevitable wear and ageing. Serious complications may result in surgical re-intervention with the possible removal of the system.

The Permedica Spa OPL system is suitable only for professional use, preferably in specialized clinical centres and in any case all user surgeons must know adequately the specific methods of application of the System.

4.2 INDICATIONS FOR USE General indications for use: ➢ Trans-femoral or trans-tibial amputation of the lower limb and trans-humeral amputation of the upper limb;

- Complications after use of external conventional socket prostheses of the lower or upper limb;
• Contraindications to the use of external conventional socket prostheses of the lower or upper limb;
• Sufficient bone quality in the stump to ensure the primary anchorage and osteointegration of the OPL intramedullary stem.

The use of devices for purposes other than those provided for is not allowed. For more detailed information, refer to the device-specific Technical Data Sheet and Surgical Technique. 4.2 CONTRAINDICATIONS Absolute contraindications include:

- Systemic infections and any septic conditions in the stump area;
• Allergy to implant component materials;
• Radiation exposure to the affected stump area rendering the bone non-viable.
Relative contraindications that could compromise the success of the intervention include:
• Chronic or acute infections of local nature, even located away from the implant site (risk of hematogenic spread of the infection to the site of the implant);
• Insufficient length of the stump to receive the implant;
• Limited articulation of the amputated limb or permanent contractures of the stump muscles;

- Vascular, muscular, neurological pathologies that put the affected anatomical part and / or the contralateral part at risk;
• Cardiovascular diseases;
• Disorders of hemostasis or clotting or healing of tissues;

- Altered bone metabolism;
• Long-term therapies of cortisone or ongoing chemotherapy treatments;
• Uncontrolled Diabetes mellitus;
• Compromised immune system such as organ transplant recipients under immunosuppression or HIV;

- Limited articulation of the amputated limb or permanent contractures of the stump muscles;
• Pregnancy, breastfeeding;
• Obesity, overweight (weight > 100 kg - for lower limb prostheses);
• Minors pre-pubertal with incomplete skeletal formation;

- Patients suffering from psychiatric disorders or mental instability that affect their compliance with the protocol;
• Patients unwilling or unable to follow the instructions for rehabilitation and aftercare provided by their doctor;
• Patients suffering from physiological instability (e.g. labyrinthitis, dizziness, etc.) which might be a risk of accidental falls;
• All concomitant pathologies that could compromise the proper functioning of the System.

4.3 RISK FACTORS The following risk factors can lead to unsatisfactory results:

- patient's inability to understand and follow the surgeon's recommendations and the physical therapy program;
• Lack of proper coupling between modular components
• Carrying out physical, work or sports activities associated with strong stresses that

belonging to the same system. Always check the compatibility of the system on the Product and / or Surgical Technical Data Sheet.

In addition, the implantation, and any removal, must be performed using the surgical instruments provided by Permedica Spa. The use of different instruments can lead to the unproper positioning of implanted components. Permedica Spa excludes any liability for the consequences deriving from the use of its devices in combination with those of other manufacturers or deriving from the use of third-party instruments. It is not permitted to reuse a device previously implanted in the patient's or another patient's body, or a device that has come into contact with another patient's body fluid or tissue, since, otherwise mechanical and / or superficial and / or geometric and / or biological integrity cannot be guaranteed (single use devices). Remove the device from the protective package respecting common aseptic techniques and manipulating it carefully. Before use, it is necessary to visually check that there are no defects. Even small scratches, dents, carvings or impurities on surfaces can lead to excessive wear and be a source of complications. Caps, capsules or any other protective equipment that may exist must be removed immediately before use. Devices with a porous titanium coating should be handled carefully to avoid damage to the coating itself. Prevent coated parts from coming into contact with materials other than the package, clean surgical gloves and patient tissues. The intramedullary stems of the OPL System must never be cemented. Avoid prolonged contact - direct or indirect - of the electrocoagulator with implantable components as it can cause structural alterations that can change the fatigue resistance characteristics of the material with consequent risk of breakage. Always remove surgical debris before closing.

## 5.5 POST-OPERATIVE STAGE

### General information for patients

The surgeon will have to inform patients in detail about the general surgical risks of the implant, the possible side effects and the measures one should take to reduce the incidence of these factors.

In particular, the patient must be informed by the surgeon regarding:

- The efforts required in terms of the length and speed of rehabilitation and the number of postoperative checks required.
- The frequency of clinical checks.
- The use and maintenance of the OPL system.
- Regular and daily cleaning of the coupling and hygiene of the stoma and skin.
- The load bearing and movements allowed during post-operative rehabilitation before the use of a lower or upper limb exoprosthesis.

Patients should be advised that any change in performance (mobility, pain, etc.) should be reported to the surgeon.

Patients must have understood all the information provided to them in order to express in a favourable and reasoned way their consent to undergo the surgery.

As a guideline Permedica Spa advises the surgeon to provide the patient with the following recommendations for the use, cleaning, hygiene and maintenance of the OPL System:

- Apply an iodopovidone-based solution daily to the skin and soft tissues around the stoma for at least 14 days after stump stoma surgery;
- Clean daily and regularly twice a day the skin and soft tissues around the stoma with drinking water and neutral soap allowing to dry in the air during the period in which the percutaneous bone fixation system is in situ and, in general, for the period in which the transcutaneous passage is present;
- Do not use disinfectants or antibiotic solutions to clean the stoma, unless directed by the surgeon;
- To remove the external prosthetic components, the patient must first release the limb exoprosthesis from the external adapter and then external coupling adapter;
- Once the exoprosthesis has been removed from the percutaneous coupling, the external adapter must be completely disassembled according to the indications provided by the orthopaedic staff or hospital treating the patient before being cleaned, unless they expressly prohibit in doing so. The individual components must be cleaned for example with hydrogen peroxide at 36 vol. or with an appropriate detergent, or with a mixture of 70% medical alcohol and 30% distilled water. After cleaning, the individual components should be rinsed in water preferably distilled and well dried with a clean cloth.

If you notice any damaged component this must in any case be replaced with a new or at least functioning component.

Patients should also be informed that any implanted metal devices:

- > could influence the result of examinations performed using ionizing fields (CT computed tomography);
- > can be detected by surveillance systems based on magnetic fields;
- > in the case of cremation, removal of the prosthesis could be required according to local provisions and regulations.

The surgeon must provide patients with the implant card filled in the required sections and advise that further information is available on the website indicated on the implant card. It is recommended to carry out regular postoperative follow-up to detect the first signs of wear, dislocation of the prosthesis, etc., and to evaluate the actions to be taken.

### MRI compatibility

The components of the OPL System have not been tested for safety and compatibility in an MR (Magnetic Resonance) environment, nor to determine their specifications for conditional status. They have not been tested for heating, migration, or image artifacts in the MR environment.

It is known that there are currently several manufacturers and different generations of MRI systems in use. Permedica Spa is not able to define the conditions under which the devices can be used in the MR environment. Therefore, the safety of implants in the MR environment is unknown, and the scan could result in injury to patients with the implant. The surgeon should make the patients aware about the risks connected to the exposure to magnetic fields. The patient must also be informed that implants can affect the results of magnetic resonance imaging (MRI) scans.

## 5.6 POSSIBLE SIDE EFFECTS

The following are the possible known side effects related to the implantation of an intramedullary percutaneous bone fixation system for external limb prostheses:

- General risks associated with anesthesia; General risks associated with major surgery;
- Hypersensitivity reactions to the metal. The significance of this sensitization needs further clinical investigation;
- Cardiovascular disorders including deep vein thrombosis, thrombophlebitis, pulmonary embolism, wound hematoma and avascular necrosis;
- Nerve damage, neuromas;
- Intraoperative bone injuries, including fissures, perforations or fractures of the diaphysis due to incorrect size, type of stem, excessive press-fit or explantation;
- Abnormal pain (chronic, persistent, psychosomatic) of the stump;
- Heterotopic ossification;
- Bacterial colonization of the stoma without signs or symptoms of infection;
- Complications in the stoma region including edema, swelling, tissue irritation, serum leakage, possible delay or failure to heal and epithelialization of the stoma;
- Superficial infections in the stoma region or percutaneous tract;
- Deep intramedullary infections, osteomyelitis, septic arthritis;
- Muscle contractures, enthesitis, tendinitis;
- Loosening of tightening between components of the prosthetic system;
- Mechanical failure or failure of the components of the prosthetic system due to excessive loads (excessive physical activity, sport, non-compliance with contraindications) or non-physiological stresses (falling to the ground and accidental trauma) or due to an incorrect choice of components (eg. undersizing of the stem and the percutaneous coupling);
- High risk of accidental fall to the ground in case of rupture of the prosthetic components;
- Possible excessive tension on the soft tissues of the stump such as skin, scars, or muscles;
- Perioperative bone fracture for postoperative trauma;
- Resorption and atrophy of the periprosthetic bone, osteonecrosis;
- Retrograde migration of stem;
- Failure or insufficient osseointegration of the intramedullary stem resulting in mobilization of the device due to e.g. poor bone quality, insufficient primary stability (press-fit), micromovements, loading and early mobilization;
- Allergic reactions to implanted materials or for metal ions released by the device;
- Complications arising from the modular connection of the prosthetic system (metallolysis, adverse biological reactions, osteolysis due to the release of particles and ions due to corrosion, abrasion, wear at the interface of the components);
- Surgical intervention to replace one or more components of the implant (e.g. percutaneous coupling before the end of the useful life, tightening screw, or the entire implant);
- Surgical intervention to permanently explant the intramedullary implant fixation;
- Surgical reintervention of stump amputation at a higher (proximal) level due to tissue necrosis, deep severe chronic infection or fracture.

## 5.7 ADVERSE EVENTS

In accordance with the definition of incident / serious incident reported in EU regulation 2017/745, any serious incident occurring in relation to the device must be reported to Permedica Spa and to the competent authority of the Member State in which the user and / or patient is established.

## 6. STERILITY

### a. Implantable medical devices

The implantable devices supplied in the "STERILE" state must remain in their closed protective packaging until the moment of implantation. Before using the device:

- check the sterility expiration date (year / month) on the product label;
- visually check that the secondary packaging and labeling are intact;
- visually check that the sterile primary packaging is intact and does not show cracks, holes or any other types of damage.

Implantable devices supplied in the "STERILE" state cannot be re-sterilized by the user. In case of accidental opening, sterility expiration date exceeded or damaged packaging, do not use the product and return it to Permedica Spa.

The sterilization method used is indicated on the label of each device. Metal implantable devices are sterilized by irradiation (minimum nominal dose of 25 kGy) or by vaporized hydrogen peroxide.

### b. Non-implantable Medical devices

Non-implantable devices supplied in the "NON-STERILE" state should not be sterilized. It is recommended to clean the product according to the recommended cleaning and maintenance procedures set out below. The performance and integrity of the product are guaranteed only if the relative packaging is not damaged. The components supplied in a "NON-STERILE" state, at the time of first use, must be properly cleaned with the normal sanitization procedures used in hospital facilities. They can be cleaned, for example, with a mixture of 70% medical alcohol and 30% distilled water. After cleaning, the product should be rinsed in sterile distilled water and well dried with a clean non-woven cloth.

### c. Surgical instruments

The surgical instruments are supplied in the "NON-STERILE" state and must be cleaned, disinfected and sterilized before use, in accordance with methods validated by the end user. For handling and sterilization of non-sterile instruments, refer to the instructions for use "Instructions for the care, cleaning, disinfection and sterilization of surgical instruments".

## 7. STORAGE, USE AND DISPOSAL CONDITIONS

All devices must be stored at room temperature and in a dry place, protected from direct exposure to sunlight or heat. They must be kept inside their original sealed protective packaging.

Once opened the package of the implantable devices, make sure that the model and size correspond to what is stated on the label. Avoid contact with objects or substances that may alter the sterile condition or surface integrity. The devices removed from the package must not be used if accidentally dropped or if they have suffered accidental trauma / impact. The code and lot number of the devices used must be recorded on the patient's medical record and on the Implant Card using the labels supplied in the package.

Disposal of implantable devices must be carried out by hospitals in accordance with applicable laws.

## 8. CUSTOM-MADE DEVICES

A custom-made device is an implantable device specifically designed and manufactured on the base of a written prescription and specific design indications provided, under their own responsibility, by any person authorized by national law by virtue of their professional qualification, in the event that, in the opinion of the medical surgeon, a standard implant cannot be used for an individual patient.

The surgeon must be aware of the limitations of the special device to be considered for the construction and choice of materials. The surgeon must also have the skills and experience that allow him to guarantee the right indications, provide the correct specifications and perform the optimal application of the customized device. Custom manufacturing may not have a corresponding instrumentation.

The custom-made devices are made according to internal procedures, based on the state of the art and the experience of Permedica Spa with standard systems, which allow to follow the development process according to the standards and regulations in force. The design and development are carried out by Permedica Spa together with the surgeon who also approves the final design proposal containing all the information of the tailored case (including data of the patient for whom it is intended). Since these are custom-made products, no laboratory or clinical tests are performed, therefore higher risks must be considered compared to a standard product.

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